



Fitting EXPERT Progressive is simple and the lenses can be ordered with K readings and Rx, or supply the lab with the information below on the EXPERT Progressive Design Order Sheet to improve fitting success with the first lens.

Please download this form, complete and fax to: **800.648.2272**

Order from:

Patient: _____ Phone: _____

Account #: _____ Email: _____

RIGHT EYE					LEFT EYE			
				Distance Rx Recent and optimized				
				Reading Add				
Flat K		Axis		K-readings	Flat K		Axis	
Steep K		Axis			Steep K		Axis	
								
				Palpebral fissure height (a)				
				Pupil diameter (b)				
				Corneal diameter (c)				
				Height from bottom of pupil to top of lower lid in mm (d)				
<ul style="list-style-type: none"> Lower lid can be at limbus (best) Lower lid can be above limbus (mark by + mm) Lower lid can be below limbus (mark by - mm; least ideal position) 				Lid position: Quantify in mm the position of the lower lid relative to the limbus 	<ul style="list-style-type: none"> Lower lid can be at limbus (best) Lower lid can be above limbus (mark by + mm) Lower lid can be below limbus (mark by - mm; least ideal position) 			
Tight	Average	Loose		Lid tonicity	Tight	Average	Loose	